

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

D 4 SURGICALS

Leading Hospital Furniture Manufacturer from India, supplying and exporting quality hospital furnitures and Surgical Equipments for over three decades.

Our Clientele

Includes many well-known hospitals ,Healthcare Institutions in India and Abroad.

Our Business

Contributing to better healthcare by designing, manufacturing, distributing and maintaining sophisticated Surgical equipments and hospital furniture is our main business.

Our Commitment

To offer better products in terms of aesthetics, reliability, simplified operation, durability, high serviceability and exceptional quality.

Our Specialisation:

- Offering total solution for small hospitals.
- Furnishing total projects, right from furniture to hospital management software.
- Specially designed projects for third world countries where cost effectiveness and functionality are prime factors in building projects .We cater to countries in Asia, Africa and central Asian republics.
- Spreading healthcare to the masses in semi-urban and small cities.
- Turning dreams to set up small hospitals into realities for doctors, medical professionals and healthcare industry companies.

Our Strength

We have a team of qualified & experienced eminent persons viz, doctors who are qualified in hospital projects and planning, architects with experience of over a decade in USA, chartered accountant having been associated with hospital projects, besides our technical personnel with varying experience in medical electronics and trade for past 2 to 4 decades.

Our mission

To provide high standards Surgical equipment and hospital furniture at attractive prices. We invite queries from the buyers and assure them a quick response.

Our Quality Assurance

We are An ISO 9001:2000 Certified Company and Member WTC Inc. New York

Our Vision

PATIENT'S COMFORT

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Electric Bed [D4S 111]

Overall Size : 2300mm X 980 mm X 475mm

Height adjustment : 475 mm to 745 mm

- Backrest movement 80 DEGREE
- Trendelenburg movement 12 DEGREE
- Collapsible railings
- Buffers at four corners
- All four movements operated by wired remote
- 150 mm dia, 4 heavy duty wheels, 2 with stopper
- **OPTIONAL: C.P.R. FACILITY, nurse control panel AND BATTERY BACKUP**



C.A.R. - Central African Republic
U.A.E. - United Arab Emirates

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

NEW CONCEPT FOR OUR MOST VALUABLE CUSTOMERS

DESIGN YOUR OWN ELECTRIC BED

STEP 1

SELECTION OF BED SURFACE



PERFORATED SHEET TOP



TRAPEZOID TYPE TOP



ABS PLASTIC TOP

STEP 2

SELECTION OF HEAD AND FOOT BOARD



ABS PLASTIC BOARDS



S.S. METAL BOWS



M.S. POWDER COATED BOWS

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

STEP 3

SELECTION OF CASTORS



DELUX CASTOR



STANDARD CASTOR

STEP 4

SELECTION OF RAILINGS



ABS PLASTIC SWING TYPE RAILING



S.S. SWING TYPE RAILING



S.S. COLLAPSIBLE RAILING



ALLUMINIUM COLLAPSIBLE RAILING



**M. S. POWDER COATED
COLLAPSIBLE RAILING**

**SELECT THE REQUIRED PARTS OF THE BED AND FORWARD US BY EMAIL, FAX
AND WE WILL QUOTE YOU THE PRICES ASAP.**

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Electric Homecare Beds

SPECIAL FEATURES :

Complete Wooden Body

- Made from Burma Teak Wood one of the highest quality of wood.
- Look like a part of your designed furniture and don't make your home look like a hospital.
- Controlled with remote control.
- Three different comfortable positions, backrest position, leg rest position and back and leg rest position.
- Size??? To suit your Bedroom.



Backrest Position



Legrest Position

U.A.E. - United Arab Emirates



Backrest + Legrest Position

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT


HOSPITAL BEDS MECHANICALLY OPERATED



Plain Hospital Bed (DELUX) [D4S 01]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Perforated CRCA sheet top
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Plain Hospital Bed (STANDARD) [D4S 01A]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Perforated CRCA sheet top
- M.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Manual Back Rest Bed [D4S 02]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Two section perforated CRCA sheet top
- Manually operated, 5 position back - rest on ratchet
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT



Semi-Fowler Bed [D4S 03]

Overall Size : L 2010 mm x W 865 mm x H 590 mm

Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Two section perforated CRCA sheet top
- Back rest operated by screw with folding handles.
- Detachable S.S. bows of 32 mm dia
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



Semi-Fowler Bed [D4S 03A]

Overall Size : L 2010 mm x W 865 mm x H 590 mm

Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Two section perforated CRCA sheet top
- Back rest operated by screw with folding handles
- Detachable S.S. bows of 32 mm dia. With metal panel (Different colour options available)
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated

CAR - Central African Republic
U.A.E. - United Arab Emirates

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Fowler Bed [D4S 04]

Overall Size : L 2060mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Four section perforated CRCA sheet top
- Back rest and knee rest operated by attached folding handles
- Detachable s.s. bows with metal panels (Different colour options available)
- S.S. swing type railings: Optional
- 125 mm dia heavy duty wheels; 2 with brakes
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Recovery Bed [D4S 05]

Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Four section perforated CRCA sheet top.
- 3 Separate screws for back rest, knee rest, trendelenburg & reverse trendelenburg position by attached folding handles.
- Detachable s.s. bows with metal panel (different colour options available).
- S.S. side railings on both sides (Swing type or collapsible type) and buffer on all corners.
- 125 mm dia heavy duty wheels; 2 with brakes.
- 4 flush I.V. rod locations.
- Pre-treated and epoxy powder coated.



D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

ICU Bed [D4S 06]

Overall Size : L 2060mm x W 905 mm x H 590 mm

Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Four section perforated CRCA sheet top
- 4 Separate screws for back rest, knee rest trendelenburg & reverse trendelenburg and height adjustment position by attached folding handles. (ALSO AVAILABLE WITH FIXED HEIGHT OPTION MODEL D4S 006-A)
- S.S. side railings on both sides (Swing type or collapsible type) and buffer on all corners
- Detachable s.s. bows with metal panel (different colour options available) . .
- 125 mm dia heavy duty wheels; 2 with brakes
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated



C.A.R. - Central Africa Republic
U.A.E. - United Arab Emirates

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT



**PLAIN HOSPITAL BED WITH CRADDLER
[D4S 001B]**

**Overall Size : L 2010 mm x W 865 mm x H 590 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube**

- Perforated CRCA sheet top
- S.S. bows of 32 mm dia.
- 4 flush I.V. rod locations
- Pre-treated and epoxy powder coated

**CRIB SPECIFICATIONS
SIZE (36" x 15") -
M.S Tubular pipe constructi
Pre-treated and Epoxy Powder coated**



PAEDIATRIC BED [D4S 000]

Size : 137cm x 76cm x 60cm.

- Complete with full-length dropside safety railings
- Head and leg bows of equal height
- Pre-treated and epoxy powder coated

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

HOSPITAL TABLES

D4S

Operation Table (Height Adjustable) [D4S 20]

Overall Size : L 1820mm x W 560 mm x Ht. Adjustable from 750mm to 920 mm

- Stainless steel top with kidney bridge
- Detachable and folding foot end
- Height Adjustment and head low/up position by screw and lever mechanism
- Pre-treated and epoxy powder coated
- Available in stainless steel top & epoxy powder coated base / S.S. base



Labour Table (Foot End Folding) [D4S 22]

Overall Size : L 1830mm x W 870mm x Ht. 800mm

- Two section Stainless steel top with head low/up position with screw mechanism
- Detachable and folding foot end
- 'U' notch with drainage tray
- Available in full stainless steel / stainless steel top & epoxy powder coated base

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Labour Table (Telescopic) [D4S 23]

Overall Size : L 1830mm x W 870mm x H 800mm

Frame : 50 mm x 25 mm Rectangular Stainless steel tube

- Three section Stainless steel top
- Back rest on ratchet; head low/up by screw mechanism
- Middle Section with 'U' notch with bucket / bowl
- Sliding foot end
- Railings on three sides
- Available in full stainless steel / stainless steel top & epoxy powder coated base
- Lithotomy Bars Gutter type / Leather Strip type



Examination cum Gynec Table [D4S 31]

Overall Size : L 1820 mm x W 560 mm x H 815 mm

- CRCA sheet / S.S. top
- 'U' notch with foot end folding and detachable
- Lithotomy bars
- Pre-treated and epoxy powder coated / Full S.S.

Patient Examination Couch [D4S 32]

Overall Size : L 1830mm x W 575 mm x H 850 mm
Frame : CRCA Sheets and Rectangular tubes

- 4 Cabinets; 1 with lock, drawer with lock
- 100MM THICK FOAM Mattress
- Sliding Stepping Stool and Writing Pad
- Pre-treated and epoxy powder coated



Gynec Examination Couch [D4S 33]

Overall Size : L 1830mm x W 575 mm x H 850 mm
Frame : CRCA Sheets and Rectangular tubes

- 3 Cabinets; 1 with lock, drawer with lock.
- Foam Mattress with 'U' notch 100 mm thick
- Examination Trolley with three trays
- Sliding Stepping Stool and Writing Pad
- Pre-treated and epoxy powder coated

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

TROLLEYS

Stretcher Trolley (Folding) [D4S 16]

Overall Size : L 1830mm - 990mm x W 570 mm x H 860 mm
Frame : 50 mm x 25 mm, CRCA Rectangular tube

- Folding Stainless steel top fixed on CRCA tube frame mounted on 150 mm dia. heavy duty wheels
- Provision for fixing I.V. rod
- Pre-treated and epoxy powder coated



Stretcher-On-Trolley [D4S 17]

Overall Size : L 1830 mm x W 530 mm x H 870 mm
Frame : 50 mm x 25 mm, CRCA tube

- Detachable Stainless steel top fixed on CRCA tube frame mounted on 150 mm dia. heavy duty wheels
- Provision for fixing I.V. rod
- Pre-treated and epoxy powder coated



D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Casualty Trolley [D4S 18]

Overall Size : L 1830 mm x W 630 mm x Ht. adjustable from 640 mm to 810mm

Frame : 50 mm x 25 mm, CRCA tube

- Height adjustable trolley with X-Ray permeable detachable top and full length X-Ray cassette holder
- Height adjustment, trendelenburg & reverse trendelenburg positions by screw and lever mechanism
- 150 mm dia. heavy duty wheels, 2 with brakes
- Provision for fixing I.V. rod
- Swing away type S.S. railings on both side, rubber buffers on all corners
- Pre-treated and epoxy powder coated



Emergency Trolley [D4S 26]

Overall Size : W 720mm x D 650mm x H 1600mm

- 6 Drawer Modular system upper two drawers with medicine containers; pre-treated and epoxy powder coated
- 6 Detachable plastic bins
- 2 cages for oxygen and nitrous cylinders
- Top tray for monitor and pulse oxymeter
- System mounted on 100 mm dia. wheels, with pair of front brakes
- Available in Full S. S. also



D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Baby Trolley [D4S 36]

Overall Size : W 390 mm x D 800 mm x H 1080 mm

- Stainless steel tray in head low position
- 3 drawer modular system
- Top with medicine containers
- Adjustable I.V. rod & lamp arrangement
- Bottom tray for mounting suction apparatus
- Trolley mounted on 100 mm. dia. wheels
- Pre-treated and epoxy powder coated



Oxygen Trolley [D4S 34]

A) Pull Type with three wheels of 50 mm dia.

B) Push type with two wheels of 100 mm dia

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT

Monitor Trolley (for scopic surgery) [D4S 29]

Overall size : W 660 mm x D 560 mm x H 1498 mm

- Available in full Stainless steel & epoxy powder coated
- Adjustable Shelves
- 125 mm dia. heavy duty castors; 2 with brakes

Food Trolley [D4S 14]

Overall Size L 760 mm x W 360 mm x Ht. adjustable from 750 mm to 1050 mm

- Stainless steel top
- Height adjustable by lever gear mechanism
- Mounted on 50 mm dia wheels
- Pre-treated and epoxy powder coated



C.A.R. - Central Africa Republic
U.A.E. - United Arab Emirates

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT



Kick Bucket and Bowl Stand [D4S 27]

- Stainless steel bowl dia - 345 mm.
- 2 Stainless steel bowls of dia - 345mm each

Stand available in stainless steel / epoxy powder coated



Instrument Trolley [D4S 24]

Overall Size : L 710mm x W 460mm x Ht. 830mm
Frame Stainless steel square tubes

- Stainless steel bottom and top shelves with railing on top shelf
- Trolley Mounted on 100 mm dia. wheels
- Available in full stainless steel / Powder coated frame with stainless steel shelves and railings



Dressing Trolley with Bowl and Bucket [D4S 35]

Overall Size : L 1080mm x B 460 mm x H 830 mm
Frame Stainless steel square tubes

- 2 stainless steel shelves with railing on top
- Stainless steel bowl dia 345 mm.
- Bucket dia. 300 mm.
- Available in full stainless steel/ Powder coated frames with stainless steel shelves and railings



Mayo's Trolley [D4S 25]

Tray size : L 555mm x W 430mm x Ht. adjustable from 860mm to 1300 mm
Frame Stainless steel rectangular and square tubes

- Trolley Mounted on 50 mm dia. wheels
- Available in full stainless steel / Powder coated frame with stainless steel tray



Multipurpose Bed Side Lockers / Cabinets [D4S 214]

Overall Size : W 410mm X D 370mm X H 900mm

- Complete ABS (Acrylonitrile Butadiene Styrene) plastic
- Different colors to select (only for bulk quantity)
- OEM supplies available by incorporating your logo (only for bulk quantity)
- Can be used in Hospitals / Hotels / Offices / Laboratories etc.



Bed Side Locker (Delux) [D4S 11]

Overall Size : W 395mm x D 305 mm x H 800 mm

- Crca sheet main cabinet and drawer
- Mounted on 50 mm dia twin roller wheels front wheels with break
- S.S. top with three sides raised and front turned down
- Pre-treated and epoxy powder coated



Bed Side Locker (standard) [D4S 10]

Overall Size : W 370mm x D 305 mm x H 800 mm

- Two Stainless steel tops with 3 sides raised and front turned down
- Railings on 3 sides above main cabinet
- Handle and magnetic latch on door
- CRCA Steel box mounted on 20 mm x 20 mm CRCA sq. tube frames
- Pre-treated and epoxy powder coated



Multipurpose Stool [D4S 28-A]

Overall Size : Top 280 mm x 280 mm x H 445 mm

- Stainless steel top / Full S.S.
- Pre-treated and epoxy powder coated



Revolving Stool [D4S 28]

Overall Size : Top 300 mm dia. with Ht. adjustable from 510 mm to 710 mm

- Stainless steel top / Full S. S.
- Pre-treated and epoxy powder coated

D 4 SURGICALS...GLOBAL VISION FOR THE PATIENT'S COMFORT



Instrument Cabinet [D4S 30]

Overall Size : W 610 mm x D 355 mm x H 1525 mm
Frame : CRCA Rectangular tube, CRCA sheet panels


- Pre-treated and epoxy powder coated /full stainless steel
- 5 Shelves, doors & 2 sides of 5 mm thick fine glass



Bed side screen [D4S 13]

Overall Size : L Open 2360 mm /
close 1200 mm x W 550 mm x H 1730 mm

- CRCA tube structure mounted on 50 mm dia wheels
- 2 Sides folding



Saline Stand [D4S 37]

Ht. adjustable from Min.:
1370 mm to Max. :2435 mm
• Available in full stainless steel /
epoxy powder coated stand & stainless
steel saline rod

Ultrasonic cleaner,
Operation theatre lamp
Wheelchairs

O. T. fumigators

Anesthesia Machine

Electrosurgical Generator

Suction Apparatus

Autoclaves

Infant Incubators

Infant Radiant Warmer

Resuscitation Unit

Phototherapy Unit

Orthopaedic Implants

Hospital Hollowwares

Surgical Instruments

Rehabilitation Products

Laboratory Equipments

Laboratory Plasticware

Laboratory Glassware

Microscopes

CY Medical Products

Herbal Cosmetics

Herbal Beauty Products

Laboratory Chemicals

**For your needs of all type of surgical and laboratory products
our sister concern RAJ INDUSTRIES is always there
to help you to find out OEM suppliers On liaison basis**

List of products we can search for you....

And many more.....

D4S

D4 SURGICALS (INDIA) PVT. LTD.

307, Samruddhi Commercial Complex Chincholi Bunder Road, Malad (West), Mumbai 400064, INDIA.

TEL : 91-22-65235614 / 15 / 40034801 / 64167621 / 28764161 / 28744280.

Fax : 91-22-28781783 / 28764161 Mobile : 9820138189 / 9821676944

email info@d4surgicals.com / exports@d4surgicals.com / sales@d4surgicals.com

WEBSITES

www.d4surgicals.com (showroom website with 12 languages)



www.surgicalport.com (different navigation style)

www.hospital-furniture-india.com (some add ones in text & products)

www.exporterhospitalfurniture.com (exclusively for dealers showing prices)

www.d4medical-surgicalequipments.COM (site on B2B portal India mart)

www.motorized-beds.com (exclusive site on motorized beds a MINI PORTAL)

d4surgicals.tradeindia.com (site on B2B portal Trade India with language converter / You can call me free)

in101415980.trustpass.alibaba.com (site on B2B portal ALIBABA.com with TRUST PASS verified)

www.exportershospitalfurniture.COM (site on B2B portal Exportersindia.com)

BRANCH OFFICE / SHOWROOM

54, OM ANAND INDUSTRIAL PREMISES, RAGHUNATH NAGAR, THANE (WEST) INDIA. PIN : 400 604

C.A.R. - Central Africa Republic
U.A.E. - United Arab Emirates

*** Due to continuous development in designs, Technical Specifications are subject to change without prior intimation.**