


Aquaject Plus

Disposable Lens Delivery System for Acrylic Foldable Lens

Contain one sterile Injector, Cartridge & Cushion



(Having Inner side Slot)



(Having Front side Slot)

Aquaject Cartridge

Contain one sterile Cartridge & Cushion



Aquaject Cartridge

Model	Recommended Incision Size
SMALL	2.2mm. to 2.8mm.
MINI	1.8mm. to 2.2mm.
Flyglide	1.7mm. to 2.40mm.
Flyglide Micro	1.4mm. to 2.0mm.



REGD. OFFICE :
OMNI LENS PVT. LTD.
5, "SAMRUDDHI", Opp. Sakar-III,
Navrangpura, Ahmedabad - 380 014. (INDIA)
Phone : +91-79-27544031-32-33-34 Fax : +91-79-27544035
E-mail : info@omnilens.in, info@bgpgroup.com
Website : www.omnilens.in



An ISO 9001 & 13485
Certified Company