

Omni - CTR

Capsular Tension Ring



Specification

Capsular Tension Ring		
Material :Clear PMMA / Blue PMMA		
Model	Size	Thickness
CTR 100	8-10 (± 0.25)	0.18 (± 0.02) mm
CTR 110	9-11 (± 0.25)	0.18 (± 0.02) mm
CTR 120	10-12 (± 0.25)	0.18 (± 0.02) mm
CTR 130	11-13 (± 0.25)	0.18 (± 0.02) mm
CTR 140	12-14 (± 0.25)	0.18 (± 0.02) mm

Description:

- One Piece PMMA ring
- High flexibility due to special design
- ETO sterilized
- MATERIAL - PMMA CQUV


Features:

- Circular expansion and stabilization of the periphery of the capsular bag
- Reduced risk of capsular fibrosis
- Reduced risk of capsular shrinking
- Safe IOL centration in eyes with zonular dehiscence

Omni

CTR INJECTOR

The Omni CTR can be comfortably inserted in the capsular bag in a very controlled manner during the surgery with the help of Omni CTR injector system, which is extremely easy to use.



REGD. OFFICE :
OMNI LENS PVT. LTD.
 5, "SAMRUDDHI", Opp. Sakar-III,
 Navrangpura, Ahmedabad - 380 014. (INDIA)
 Phone : +91-79-27544031-32-33-34 Fax : +91-79-27544035
 E-mail : info@omnilens.in, info@bgpgroup.com
 Website : www.omnilens.in



An ISO 9001 & 13485
 Certified Company