


CE
1293



The Ideal Choice
For Precision Control & Consistency In Small Incision IOL Implantation




MICRO SURGICAL KNIVES FOR SMALL INCISION SURGERY




◀ **LANCE TIP KNIFE** – For first scleral Incision and can also be used for making side port

- 15° Blade, Straight




◀ **CRESCENT KNIFE** – For making groove or tunnel incision

- Straight / Angled Bevel up/down




◀ **SLIT KNIFE (POINTED TIP)** – For entry in to Anterior Chamber

- 2.5mm Straight/Angled Bevel up/down
- 2.8mm Straight/Angled Bevel up/down
- 3.0mm Straight/Angled Bevel up/down
- 3.2mm Straight/Angled Bevel up/down



◀ **SLIT KNIFE (BLUNT TIP)** – For extension of Slit.

- 5.1mm Straight/Angled Bevel up/down
- 5.2mm Straight/Angled Bevel up/down
- 5.5mm Straight/Angled Bevel up/down



◀ **MVR KNIFE** – For side port incision and vitreoretinal Procedures.

- 19 G round Shaft, double cutting edge
- 20 G round Shaft, double cutting edge
- 23 G round Shaft, double cutting edge
- 25 G round Shaft, double cutting edge

Features :

- Available in Flat & Round Stock.
- Manufactured from High Quality Steel.
- Ready to use, ETO Sterile.
- Made from Poly Carbonate material.
- Manufactured under high magnification.
- Custom-design facility available.



REGD. OFFICE :
OMNI LENS PVT. LTD.
5, "SAMRUDDHI", Opp. Sakar-III,
Navrangpura, Ahmedabad - 380 014. (INDIA)
Phone : +91-79-27544031-32-33-34 Fax : +91-79-27544035
E-mail : info@omnilens.in, info@bgpgroup.com
Website : www.omnilens.in

An ISO 13485 Certified Company



CE
1293