

Omni Visc (Pfs / Vial)

Hydroxypropyl Methyl Cellulose Ophthalmic Solution Usp



(Supplied 2 ml. in a sterile pre-filled glass Syringe with grip and a sterile single use Cannula.)



5 ml.



3 ml.

(Also available in 3ml / 5ml. Vials.)

Description :

Omni Visc is a sterile, isotonic, nonpyrogenic viscoelastic solution of highly purified, non-inflammatory, 2% hydroxypropyl methyl cellulose (HPMC) with a high molecular weight greater than 80,000 daltons.

Characteristics :

Omni Visc is an ophthalmic surgical aid for use in anterior segment surgery. Due to viscoelastic properties, Omni Visc is to protect the corneal endothelium from possible damage. Used as a clinical aid to keep ocular tissues separated. Easily removable after surgery, thus minimizing the risk of trabecular obstruction that could increase the intraocular pressure following surgery.

Indications :

1. Intra Ocular Lens implant
2. Keratoplasty
3. Glaucoma Surgery

How supplied :

Omni Visc supplied in 2ml.prefilled glass syringe or 3ml./5ml in vial.

Storage : Room temperature, no refrigeration is required.



REGD. OFFICE :
OMNI LENS PVT. LTD.
5, "SAMRUDDHI", Opp. Sakar-III,
Navrangpura, Ahmedabad - 380 014. (INDIA)
Phone : +91-79-27544031-32-33-34 Fax : +91-79-27544035
E-mail : info@omnilens.in, info@bgpgroup.com
Website : www.omnilens.in

An ISO 13485 Certified Company

